

## Instructions for use of Nosten2®

### 1/ Decontaminate Nosten2®

#### 2/ Inflate the Soufflet®

#### 3/ Adjust the pressure

### 1/ Decontaminate Nosten2®

- > Put the sliding mass on the smallest number.
- > Raise and tilt the mobile pan.
- > Soak a disposable wipe with the decontaminant commonly used in the Service (e.g: Anios, Surfianos, ...).
- > Rub all the surfaces of **Nosten2®** for a perfect decontamination.

### 2/ Inflate the Soufflet®

- > Take the **Soufflet®** out of its packaging and place it between the 3 **Nosten2®** studs.
- > Take a 60cc syringe with “Luer lock” tip.
- > Screw the syringe on the “3-way stopcock” of the **Soufflet®**.
- > Inflate the **Soufflet®** by turning the tap to search for air without unscrewing the syringe.
- > Repeat the action 5 times to inject 300cc of air. The top of the **Soufflet®** should be spherical but flexible under finger pressure.
- > Replace the **Nosten2®** mobile pan on the **Soufflet®**.

The base of the front edge of the mobile pan must be at the height of the notches of Nosten2®. If it is not aligned, remove air from the Soufflet® with the syringe still connected to the “3-way stopcock”.



### 3/ Adjust the pressure

Adjust the pressure of the **Soufflet®** by positioning the sliding mass of **Nosten2®** on one of the values indicated on the top of the pan. These values correspond to 1 cm of water to the values given by the manometer.

- ⇒ For adult patients: place the sliding mass on the value 24 (*front part at the level of the notch of the sliding mass*).
- ⇒ In pediatrics: leave the **Nosten2®** sliding mass in position on the minimum pressure (*front part at the level of the notch of the sliding mass*).

Connect the “Luer lock” located at the end of the **Soufflet®** tubing to the non-return valve of the patient’s balloon air line. Adjust the pressure according to the patient’s needs by moving the **Nosten2®** sliding mass.

The correct position of the mobile pan must be checked once a day by the nursing team. If some additional air has been necessary to maintain the correct position, the other nursing teams must be informed.

#### Maintenance of Nosten2®:

To facilitate the movement of the sliding mass if necessary: lift the mobile pan, hold it vertically and place silicone oil at the top of the groove traced in the thickness of the pan.